



## **KCARE ECONOMY SHOWER CHAIR**

\$110.00 Original price was: \$110.00.<u>\$99.00</u>Current price is: \$99.00.



**SKU:** 12028S

## **PRODUCT DESCRIPTION**

# **PURPOSE & BENEFITS**

The **KCare Economy Shower Chair** provides a safe and stable seating option for users who require additional support while showering. Its strong steel frame and supportive armrests assist with positioning and transfers, helping to reduce the risk of slips and falls in the bathroom.

### **FEATURES**

- Provides secure seating while showering
- Armrests offer extra stability and assistance with sitting and standing
- Constructed from durable 22mm powder-coated white steel tubing
- Solid plastic seat for strength and easy cleaning
- Height adjustable legs for personalised comfort
- Designed and tested to Australian standards for safety and reliability

#### **SPECIFICATIONS**

	Specification	Detail
Min Overall Height	610 mm	
Max Overall Height	730 mm	
Seat Width	460 mm	



Specification	Detail
Seat Depth	480 mm
Min Seat Height	430 mm
Max Seat Height	550 mm
Maximum User Weight	110 kg