



KCARE ELECTRIC BATH LIFT



SKU: KA610

PRODUCT DESCRIPTION

Purpose & Benefits of KCare Electric Bath Lift

The KCare Electric Bath Lift is designed to provide a safe and convenient bathing experience for individuals with mobility limitations or those who require assistance getting in and out of the bathtub. This electric bath lift offers several features that promote independence, comfort, and safety during bathing.

Features

- **Powered Lifting Mechanism:** Utilises an electric motor to smoothly lower and raise the user into and out of the bathtub, eliminating the need for manual lifting or assistance from caregivers.
- **Safe and Stable Design:** Features a stable base and sturdy construction to ensure the user's safety and stability during transfers in and out of the bath.
- **Comfortable Seat:** Equipped with a comfortable and ergonomic seat that provides support and comfort while bathing, reducing strain on the user's body.
- **Waterproof Construction:** Made from waterproof materials that are resistant to moisture and corrosion, ensuring durability and longevity in wet environments.
- **Easy-to-Use Controls:** Includes simple-to-operate controls that allow users to adjust the height and position of the bath lift with ease, providing personalised comfort and support.
- **Compact and Portable:** Designed to be lightweight and portable, making it easy to transport and store when not in use.
- **Battery Operation:** Powered by a rechargeable battery, offering flexibility in placement and ensuring uninterrupted operation during bath time.
- **Safety Features:** Equipped with safety features such as a low battery indicator, emergency

stop button, and safety lockout function to prevent accidental operation.

Specifications

- **Safe Working Load:** 140kg
- **Dimensions:**
 - Overall Length 728 mm
 - Overall Width 360 mm
 - Min Overall Height 650 mm
 - Max Overall Height 1071 mm
 - Seat Width 350 mm
 - Seat Depth 430 mm
 - Min Seat Height 60 mm
 - Max Seat Height 480 mm
 - Folded Width 360 mm
 - Backrest Height 645 mm
 - Product Weight 11 kg
- **Maintenance:** Simple maintenance requirements, including regular cleaning and inspection to ensure proper function and hygiene.

The KCare Electric Bath Lift offers a safe, comfortable, and convenient bathing solution for individuals with mobility challenges, allowing them to enjoy a relaxing and independent bathing experience. With its user-friendly design and advanced features, this bath lift enhances safety and accessibility in the bathroom, promoting overall well-being and quality of life for users and caregivers alike.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au