



KCARE ELECTRIC TILT IN SPACE MOBILE SHOWER COMMODE

\$7,196.00 - \$7,785.00Price range: \$7,196.00 through \$7,785.00



SKU: KA1150

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The KCare Electric Tilt In Space Mobile Shower Commode is an innovative solution for those with mobility impairments who require assistance with showering and toileting. With electric tilt and recline functions, this shower commode is designed for maximum comfort and convenience.

The chair features swingaway footrests, height adjustable armrests, and an adjustable headrest to accommodate a variety of users. The padded seat and backrest provide added comfort and support, while the four total locking castors ensure stability and safety during use.

The electric tilt function allows the chair to be adjusted to a 30-degree angle for easy access and comfort during transfers. The chair can also be reclined up to 30 degrees for added comfort and convenience during showering or toileting.

The KCare Electric Tilt In Space Mobile Shower Commode is easy to clean and maintain, with a corrosion-resistant frame and removable, easy-to-clean seat and backrest.

Please note: Padded seat is not included

FEATURES

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- Dual electric actuator tilt mechanism.
- Contoured adjustable clamp on the headrest.
- Curved vinyl-covered padded backrest.
- Integrated push-handle for manoeuvrability.



- Curved removable padded armrests.
- Swingaway and removable height adjustable footrest with contoured padded calf support.
- 3x 125mm total locking twin castors 1x 125mm directional locking twin castor.
- Designed to be used with K Care patented removable seats.
- Available in 600 and 650mm seat widths.
- Seat, pan, and pan carrier not included sold separately

SPECIFICATIONS

Overall Width 755 mm
Overall Depth 1140 mm
Overall Height 1230 mm
Seat Depth 480 mm
Floor To Seat Height 550 mm

Armrest Type Padded Removable
Armrest Height 315 mm
Backrest Height 670 mm
Footrest Type Swing Away

Product Weight 38 kg
Max User Weight 300 kg