



\*Seat not included

## KCARE ELECTRIC TILT IN SPACE MOBILE SHOWER COMMODE



**SKU:** KA1150

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The KCare Electric Tilt In Space Mobile Shower Commode is an innovative solution for those with mobility impairments who require assistance with showering and toileting. With electric tilt and recline functions, this shower commode is designed for maximum comfort and convenience.

The chair features swingaway footrests, height adjustable armrests, and an adjustable headrest to accommodate a variety of users. The padded seat and backrest provide added comfort and support, while the four total locking castors ensure stability and safety during use.

The electric tilt function allows the chair to be adjusted to a 30-degree angle for easy access and comfort during transfers. The chair can also be reclined up to 30 degrees for added comfort and convenience during showering or toileting.

The KCare Electric Tilt In Space Mobile Shower Commode is easy to clean and maintain, with a corrosion-resistant frame and removable, easy-to-clean seat and backrest.

**Please note: Padded seat is not included**

### FEATURES

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- Dual electric actuator tilt mechanism.
- Contoured adjustable clamp on the headrest.
- Curved vinyl-covered padded backrest.
- Integrated push-handle for manoeuvrability.

- Curved removable padded armrests.
- Swingaway and removable height adjustable footrest with contoured padded calf support.
- 3x 125mm total locking twin castors 1x 125mm directional locking twin castor.
- Designed to be used with K Care patented removable seats.
- Available in 600 and 650mm seat widths.
- **Seat, pan, and pan carrier not included - sold separately**

## SPECIFICATIONS

Overall Width	755 mm
Overall Depth	1140 mm
Overall Height	1230 mm
Seat Depth	480 mm
Floor To Seat Height	550 mm
Armrest Type	Padded Removable
Armrest Height	315 mm
Backrest Height	670 mm
Footrest Type	Swing Away
Product Weight	38 kg
Max User Weight	300 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)