

KCARE ADJUSTABLE HEIGHT BED CRADLE

\$88.00



SKU: KA545



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The KCare adjustable height bed cradle is designed to fit under the mattress to form a support to keep bedclothes clear of the feet and legs. The cradle is constructed so that one side is placed underneath the mattress and the other positioned horizontally above the bed holding the weight of the bed covers off the body.

Height adjustable from 480 to 700 mm, and made from powder coated steel, this bed cradle ensures you can adjust to suit both the bed and users' needs.

FEATURES

- Quick and easy to fit or remove
- Fits under the mattress
- Keeps bedding & clothes clear of feet and legs

SPECIFICATIONS

Overall Width	510 mm
Overall Length	310 mm
Min Overall Height	480 mm
Max Overall Height	700 mm