

KCARE ADJUSTABLE HEIGHT BED CRADLE



SKU: KA545

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The KCare adjustable height bed cradle is designed to fit under the mattress to form a support to keep bedclothes clear of the feet and legs. The cradle is constructed so that one side is placed underneath the mattress and the other positioned horizontally above the bed holding the weight of the bed covers off the body.

Height adjustable from 480 to 700 mm, and made from powder coated steel, this bed cradle ensures you can adjust to suit both the bed and users' needs.

FEATURES

- Quick and easy to fit or remove
- Fits under the mattress
- Keeps bedding & clothes clear of feet and legs

SPECIFICATIONS

Overall Width	510 mm
Overall Length	310 mm
Min Overall Height	480 mm
Max Overall Height	700 mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au