



## KCARE FLINDERS CARBON FIBRE WALKER



**SKU:** KA380CFS-H

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **KCare Flinders Carbon Fibre Walker** is a lightweight and durable mobility aid designed for everyday use, providing stability and comfort for individuals who require support while walking. Its carbon fibre construction offers strength while keeping the frame exceptionally light, making it easier to manoeuvre both indoors and outdoors.

Available in small, medium, and large sizes, the walker accommodates a range of user heights. The ergonomic handgrips, adjustable handle height, and supportive seat with backrest enhance user comfort. Large wheels ensure smoother travel across varied surfaces, while the foldable design allows convenient storage and transport.

### Features

- Lightweight carbon fibre frame for easy handling
- Ergonomic handgrips for user comfort
- Adjustable handle height for customised fit
- Padded seat and backrest for resting support
- Large wheels for smoother movement on various surfaces
- Foldable design for easy storage and transport

## Specifications

- **Frame Material:** Carbon fibre
- **Wheel Size:** 10 inches
- **Seat Width:** 460 mm
- **Seat Height:** Small 485 mm | Medium 560 mm | Large 610 mm
- **Handle Height:** Adjustable
- **Weight:** 5.3 kg
- **Maximum User Weight:** 135 kg
- **Manufacturer:** KCare
- **Model:** Flinders

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)