


**KCARE FOLDING WALKING FRAME - HEIGHT ADJUSTABLE****\$121.00****SKU: KA353W****PRODUCT DESCRIPTION****PURPOSE & BENEFITS**

A lightweight, height-adjustable, folding walking frame. Made from anodised aluminium tubing, with 4 large rubber ferrules & soft foam rubber handgrips.

**FEATURES**

- Lightweight aluminium walking frame.
- Height adjustable legs.
- Simple easy to use folding mechanism.
- Front wheels and rear skis allow the frame to be pushed rather than lifted.
- 5" front wheels with slides on rear.

**SPECIFICATIONS**

Overall Width	625 mm
Min Overall Height	795 mm
Max Overall Height	970 mm
Product Weight	5 kg
Max User Weight	136 kg