



# KCARE FOLDING WALKING FRAME - HEIGHT ADJUSTABLE

\$121.00 Original price was: \$121.00.\$109.00 Current price is: \$109.00.



**SKU:** KA353W

## **PRODUCT DESCRIPTION**

#### **PURPOSE & BENEFITS**

A lightweight, height-adjustable, folding walking frame. Made from anodised aluminium tubing, with 4 large rubber ferrules & soft foam rubber handgrips.

## **FEATURES**

- Lightweight aluminium walking frame.
- Height adjustable legs.
- Simple easy to use folding mechanism.
- Front wheels and rear skis allow the frame to be pushed rather than lifted.
- 5" front wheels with slides on rear.

#### **SPECIFICATIONS**

Overall Width	625 mm
Min Overall Height	795 mm
Max Overall Height	970 mm
Product Weight	5 kg
Max User Weight	136 kg