



## KCARE FOREARM CRUTCH WITH ERGONOMIC HANDLE

**\$59.00**



**SKU: KA344**

## PRODUCT DESCRIPTION

### Purpose & Benefits of the KCare Forearm Crutch with Ergonomic Handle

The KCare Forearm Crutch with Ergonomic Handle is designed to provide comfortable and reliable support for individuals requiring assistance with mobility. Its ergonomic design and adjustable features ensure optimal user comfort and adaptability.

### Features

- **Lightweight Aluminium Frame** – Constructed from lightweight aluminium, making the crutch easy to handle and manoeuvre.
- **Ergonomic Handle** – Designed to add comfort for users, reducing hand and wrist strain during use.
- **Independently Adjustable Cuff and Shaft Heights** – Allows customisation to suit each user's specific needs and preferences.
- **Non-Slip Crutch Tips** – Equipped with non-slip tips to enhance safety by providing better traction on various surfaces.

### Specifications

- **Maximum Handle Height:** 850 mm
- **Minimum Handle Height:** 620 mm
- **Maximum Overall Height:** 1350 mm
- **Minimum Overall Height:** 890 mm
- **Maximum User Weight:** 136 kg

The KCare Forearm Crutch with Ergonomic Handle offers a practical and comfortable solution for

individuals seeking enhanced mobility support.