



# KCARE FOREARM CRUTCH WITH ERGONOMIC HANDLE

\$59.00 Original price was: \$59.00.<u>\$53.00</u>Current price

is: \$53.00.



**SKU:** KA344

### **PRODUCT DESCRIPTION**

## Purpose & Benefits of the KCare Forearm Crutch with Ergonomic Handle

The KCare Forearm Crutch with Ergonomic Handle is designed to provide comfortable and reliable support for individuals requiring assistance with mobility. Its ergonomic design and adjustable features ensure optimal user comfort and adaptability.

#### **Features**

- **Lightweight Aluminium Frame** Constructed from lightweight aluminium, making the crutch easy to handle and manoeuvre.
- **Ergonomic Handle** Designed to add comfort for users, reducing hand and wrist strain during use.
- Independently Adjustable Cuff and Shaft Heights Allows customisation to suit each user's specific needs and preferences.
- Non-Slip Crutch Tips Equipped with non-slip tips to enhance safety by providing better traction on various surfaces.

## Specifications

Maximum Handle Height: 850 mm
Minimum Handle Height: 620 mm
Maximum Overall Height: 1350 mm
Minimum Overall Height: 890 mm
Maximum User Weight: 136 kg

The KCare Forearm Crutch with Ergonomic Handle offers a practical and comfortable solution for



individuals seeking enhanced mobility support.