




KCARE FUSION PRESSURE CARE CHAIR



SKU: 18030

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **KCare Fusion High Back Pressure Chair** is designed to provide exceptional comfort, posture support, and pressure relief for individuals who spend extended periods seated. Its contoured seat, high backrest, and memory foam head pillow ensure ergonomic support, while the four-way stretch fabric enhances pressure distribution and helps prevent skin breakdown. Ideal for healthcare, aged care, or home use, it promotes comfort, durability, and hygiene.

FEATURES

- Contoured seat and high backrest provide superior comfort and support
- Removable memory foam head pillow for personalised positioning
- Height-adjustable seat for user customisation
- Fusion four-way stretch fabric is waterproof, anti-bacterial, and shear-resistant
- Foam construction helps reduce pressure points and improve circulation
- Easy-care materials designed for long-term use and easy cleaning
- Suitable for healthcare, aged care, and home environments

SPECIFICATIONS

- **Seat Width:** 520 mm
- **Seat Depth Range:** 500 – 650 mm
- **Seat Height Range:** 450 – 550 mm
- **Backrest Height:** 1130 mm
- **Max User Weight:** 250 kg
- **Frame Type:** Height adjustable
- **Upholstery:** Fusion four-way stretch, waterproof, and anti-bacterial fabric
- **Padding:** Foam with removable memory foam head pillow

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au