



KCARE FUSION PRESSURE CARE CHAIR

\$1,472.00 Original price was:

\$1,472.00.<u>\$1,324.80</u>Current price is: \$1,324.80.



SKU: 18030

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **KCare Fusion High Back Pressure Chair** is designed to provide exceptional comfort, posture support, and pressure relief for individuals who spend extended periods seated. Its contoured seat, high backrest, and memory foam head pillow ensure ergonomic support, while the four-way stretch fabric enhances pressure distribution and helps prevent skin breakdown. Ideal for healthcare, aged care, or home use, it promotes comfort, durability, and hygiene.

FEATURES

- Contoured seat and high backrest provide superior comfort and support
- Removable memory foam head pillow for personalised positioning
- Height-adjustable seat for user customisation
- Fusion four-way stretch fabric is waterproof, anti-bacterial, and shear-resistant
- Foam construction helps reduce pressure points and improve circulation
- Easy-care materials designed for long-term use and easy cleaning
- · Suitable for healthcare, aged care, and home environments



SPECIFICATIONS

• Seat Width: 520 mm

Seat Depth Range: 500 - 650 mm
Seat Height Range: 450 - 550 mm

Backrest Height: 1130 mmMax User Weight: 250 kg

• Frame Type: Height adjustable

• **Upholstery:** Fusion four-way stretch, waterproof, and anti-bacterial fabric

• Padding: Foam with removable memory foam head pillow