



KCARE HANDMADE CLOSED SEAT



SKU: KA11HCCP-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **KCare Handmade Closed Seat in Blue Vinyl** is designed for superior hygiene, comfort, and durability in clinical and home care environments. It features a waterproof, antibacterial, and antifungal vinyl cover that's easy to clean and maintain. The high-quality craftsmanship ensures long-lasting use, making it ideal for commodes and shower chairs. This seat offers a secure, comfortable surface with excellent tear strength for enhanced safety.

FEATURES

- Single-density vinyl covered seat with smooth, durable finish
- Waterproof cover for easy cleaning and infection control
- Antifungal, antibacterial, and flame-retardant material
- Tested for abrasion resistance and superior tear strength
- Comfortable and hygienic for long-term use
- Available in a range of widths to fit multiple frame sizes
- Handmade for consistent quality and precision

SPECIFICATIONS

- **Overall Width:** 445 mm
- **Overall Depth:** 435 mm
- **Overall Height:** 55 mm
- **Aperture Width:** 200 mm
- **Product Weight:** 2 kg
- **Brand:** KCare

Available Sizes

- Handmade Closed Seat Blue Vinyl – 320 mm Wide
- Handmade Closed Seat Blue Vinyl – 400 mm Wide
- Handmade Closed Seat Blue Vinyl – 445 mm Wide
- Handmade Closed Seat Blue Vinyl – 500 mm Wide
- Handmade Closed Seat Blue Vinyl – 550 mm Wide
- Handmade Closed Seat Blue Vinyl – 600 mm Wide
- Handmade Closed Seat Blue Vinyl – 650 mm Wide
- Handmade No Aperture Seat Blue Vinyl – 445 mm Wide

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au