



# KCARE HEAVY DUTY SHOWER CHAIR - HEIGHT ADJUSTABLE



SKU: KA220AAHD

## **PRODUCT DESCRIPTION**

## **PURPOSE AND BENEFITS**

Designed for use in the shower, this KCare Shower Chair has a seat area of 48cm (w) x 43cm (d).

It has a one-piece contoured plastic back rest / seat to provide positive postural support and comfort. The seat is removable for cleaning and replacement. The chair is angle and height adjustable. There are arms to provide support to users while sitting or transferring. Rubber tips on each leg provide stability.

## **FEATURES**

- Heavy duty shower chair
- 1 piece contoured plastic back rest / seat
- Seat is easily removed for cleaning
- Seat height & angle can be adjusted
- Rubber tips on legs to provide anti-slip stability
- Maximum weight capacity of 160kg

## **SPECIFICATIONS**

Code KA220AAHD

Name Aluminium Heavy Duty Shower Chair

Overall Width 530mm

Overall Depth 520mm



Seat Height 430 - 550mm

Seat Width 440mm

Seat Depth 465mm

Width Between Arms 490mm

Product Weight 3.8kg

Max User Weight 200kg

For further information go to medimart.com.au or give us a call: 1300 633 063

## SA

## **Brighton Showroom**

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

#### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

## **Reynella Showroom**

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

## **VIC**

#### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

#### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au

## NT

## **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au