



KCARE HEAVY DUTY HEIGHT ADJUSTABLE SHOWER STOOL

\$216.00



SKU: KA222AAHD

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Designed to assist those who have difficulty showering. Height adjustable legs. Plastic seat for ease of cleaning. Powder coated white 22mm steel tubing.

FEATURES

- Adjustable height and angle to suit the individual user and uneven floors
- Moulded plastic clip-on seat
- Built in armrests providing user support when sitting or transferring
- Anti-slip rubber tips on each leg for stability
- Durable and long lasting with a zinc treated mild steel frame and impervious plastic seat

SPECIFICATIONS

Max Seat Height	640mm
Max User Weight	160kg
Min Seat Height	515mm
Overall Depth	380mm
Overall Width	565mm
Product Weight	4kg
Seat Depth	290mm
Seat Width	340mm