



KCARE HIGH LOW SINGLE BABY BATH

\$7,347.00



SKU: FN20-3860

PRODUCT DESCRIPTION

Purpose & Benefits of the K Care High Low Single Baby Bath

The K Care High Low Single Baby Bath is designed to provide a safe and ergonomic bathing solution for infants, ensuring comfort for both the baby and the caregiver. Its electric height adjustment feature allows caregivers to set the bath at an optimal working height, reducing strain and enhancing safety during bathing routines. This bath is ideal for use in clinical settings, nurseries, and home environments where infant care is a priority.

Features

- **Electric Height Adjustment:** Enables effortless raising and lowering of the bath to accommodate caregiver comfort and accessibility.
- **Foot Pedal Operation:** Allows hands-free height adjustments, maintaining hygiene and ease of use.
- **Wall-Mountable Battery Charger:** Ensures the bath is always ready for use with a conveniently located charging system.
- **Drain Plug, Hose, and Valve:** Facilitates easy and efficient drainage after bathing, simplifying the cleaning process.
- **IV Pole Socket:** Provides the option to attach an IV pole (not included), accommodating medical needs during bathing.
- **100 mm Diameter Castors with Locking Mechanism:** Offers mobility for positioning the bath as required, with locking castors to secure it in place during use.
- **Powder-Coated Finish:** Ensures durability and ease of cleaning, maintaining a hygienic environment.

Specifications

Overall Width	925 mm
Overall Depth	405 mm
Min Overall Height	950 mm
Max Overall Height	1310 mm
Max User Weight	100 kg