

KCARE HILITE HIGH BACK CHAIR

\$399.38



SKU: KA582KDG05



PRODUCT DESCRIPTION

Purpose & Benefits of KCare Hilite High Back Chair

The **KCare Hilite Chair** is a high-back day chair designed to provide superior comfort, postural support, and ease of use for individuals in home care, hospital, or aged care settings. Its supportive design makes it ideal for users who spend extended periods seated and require added stability and adjustability.

This knock-down version is easy to assemble and disassemble, making it convenient for transport, storage, or relocation. The chair features a contoured backrest, padded armrests, and an easy-clean vinyl finish, all in a modern greystone colour that complements any environment.

Key Features:

- High-back design for optimal comfort and postural support
- Knock-down frame for easy transport and storage
- Padded armrests for enhanced comfort

- Durable, easy-clean vinyl upholstery in contemporary greystone
- Adjustable leg height to suit individual user needs
- Non-slip rubber tips for added stability
- Ideal for home care, hospitals, and aged care facilities

Specifications:

Min Overall Width	640 mm
Max Overall Width	700 mm
Min Seat Width	490 mm
Max Seat Width	550 mm
Min Seat Depth	480 mm
Max Seat Depth	500 mm
Min Seat Height	430 mm
Max Seat Height	580 mm
Backrest Height	700 mm
Max User Weight	200 kg
Brand	K Care