



ALIVIO KAHLO LIFT RECLINER - STONEWASH FUDGE



SKU: KA557

PRODUCT DESCRIPTION

Purpose & Benefits

The **Kahlo Lift Recliner** seamlessly combines comfort, mobility, and support in a modern design. Its three-motor system allows independent adjustment of the headrest, backrest, and footrest, enabling users to find the perfect sitting or resting posture. The dual-stage vertical lift function offers smooth, safe transitions from sit to stand, especially useful for individuals with mobility challenges. With a lay-flat recline capability and footrest extension, it also delivers pressure relief and sleeping comfort. The durable Stonewash Fudge upholstery is soft and resilient, and the hand controller with preset and all-lift/all-close buttons ensures effortless positioning for everyday ease and wellbeing.

Please note: May require assembly.

Features

- Three independent motors for precise control of headrest, backrest, and footrest
- Dual-stage vertical lift assists safe and upright stand assistance
- Lay-flat recline with footrest extension for full-length pressure relief and rest
- Soft-grip hand controller with multiple preset positions, plus all-lift and all-close functions
- Stonewash Fudge fabric: ultra-durable, stylish, and easy to keep clean

Specifications

- **Overall Dimensions (L x W x H):** Approx. 183 cm x 89 cm x 109 cm
- **Seat Size (W x D):** Approx. 52 cm x 52 cm
- **Seat Height (floor to seat):** Approx. 50.8 cm
- **Safe Working Load:** Up to 158 kg
- **Product Weight:** Approx. 74 kg

[KCare Alivio User Manual](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au