


**ALIVIO KAHLO LIFT RECLINER - STONEWASH FUDGE****\$2,617.00****SKU: KA557****PRODUCT DESCRIPTION****Purpose & Benefits**

The **Kahlo Lift Recliner** seamlessly combines comfort, mobility, and support in a modern design. Its three-motor system allows independent adjustment of the headrest, backrest, and footrest, enabling users to find the perfect sitting or resting posture. The dual-stage vertical lift function offers smooth, safe transitions from sit to stand, especially useful for individuals with mobility challenges. With a lay-flat recline capability and footrest extension, it also delivers pressure relief and sleeping comfort. The durable Stonewash Fudge upholstery is soft and resilient, and the hand controller with preset and all-lift/all-close buttons ensures effortless positioning for everyday ease and wellbeing.

**Features**

- Three independent motors for precise control of headrest, backrest, and footrest
- Dual-stage vertical lift assists safe and upright stand assistance
- Lay-flat recline with footrest extension for full-length pressure relief and rest
- Soft-grip hand controller with multiple preset positions, plus all-lift and all-close functions
- Stonewash Fudge fabric: ultra-durable, stylish, and easy to keep clean

## Specifications

- **Overall Dimensions (L × W × H):** Approx. 183 cm × 89 cm × 109 cm
- **Seat Size (W × D):** Approx. 52 cm × 52 cm
- **Seat Height (floor to seat):** Approx. 50.8 cm
- **Safe Working Load:** Up to 158 kg
- **Product Weight:** Approx. 74 kg