



KCARE KARMA FOLDABLE HEADREST

\$376.00



SKU: 2-KHRF-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **KCare Karma Foldable Headrest** provides essential head and neck support for wheelchair users, folding conveniently with the chair to preserve functionality and ease of use. Height, depth, and angle adjustments ensure customised positioning, while its compact aluminium frame and breathable padding deliver durable comfort—making it ideal for daily mobility.

Features

- Folds automatically with wheelchair—no need to remove
- Height, depth, and angle adjustable for personalised head support
- 360° rotation on a ball joint (Medium variant) for complete positioning flexibility
- Lightweight aluminium frame with breathable padding ensures comfort and longevity
- Universal compatibility with most Karma manual wheelchairs

Variations & Specifications

- **Small Size**
 - **Seat Width Compatibility:** 36 cm – 46 cm
- **Large Size**

- **Seat Width Compatibility:** 46 cm – 56 cm