



KCARE MAXI OVER TOILET FRAME - HEIGHT ADJUSTABLE

\$520.00



SKU: KA410Z60

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The K Care Maxi Over Toilet Aid is designed to help larger individuals who have trouble using a standard toilet. It features a height & angle adjustable lift-up seat with splash guard, which are both removable for easy cleaning, rubber tips on each leg for maximum stability, and sturdy arms to provide extra support to the user while in use or during transfers.

FEATURES

- Reinforced frame for maximum safety
- Adjustable height and angle to suit the occupant and uneven floors
- Arms to provide support when sitting or transferring
- Removable lift-up plastic toilet seat & splash guard for easy cleaning
- Moulded PVC seat
- Rubber tips on each leg provides anti-slip stability
- Epoxy powder coated zinc treated mild steel frame

SPECIFICATIONS

Seat Width	600 mm
Seat Depth	410 mm
Min Seat Height	440 mm
Max Seat Height	560 mm
Max User Weight	250 kg

