



KCARE MAXI SHOWER CHAIR

\$800.00



SKU: KA220ZA65-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Designed for bariatric occupants to use in the shower, these KCare Maxi Shower Chairs have a higher weight capacity of 300 Kg and a larger seat area of 550/650 W x 470 D mm.

They have a padded backrest and padded seat to provide support and comfort. There are arms with padded armrests to provide support to users while sitting or transferring. Rubber tips on each leg provide anti-slip stability.

The shower chair is durable with a reinforced heavy-duty epoxy powder coated zinc treated mild steel, white frame with waterproof commercial grade vinyl upholstery.

FEATURES

- Reinforced heavy-duty frame for improved safety
- Padded backrest for extra comfort
- Arms to provide support when sitting or transferring
- Padded armrests for extra comfort
- Padded seat for extra comfort
- Larger seat area
- Epoxy coated zinc treated mild steel frame

SPECIFICATIONS

KA220ZA55

Maximum Occupant Weight 300 Kg

KA220ZA65

Width	685 mm	785 mm
Depth	625 mm	
Height	940 mm	
Width Between Arms	565 mm	665 mm
Height of Arms from Seat	170 mm	
Back Rest Height	420 mm	
Seat Width	550 mm	650 mm
Seat Depth	470 mm	
Seat Height	515 mm	
FINISH		
Frame	Epoxy powder coated zinc treated mild steel, white	
Back Rest	Commercial grade vinyl upholstery, white	
Seat	Commercial grade vinyl upholstery, white	