



KCARE TRANSIT MOBILE SHOWER COMMODE - HEIGHT ADJUSTABLE

\$1,598.00



SKU: KA114S-HA

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The KCare Mobile Height Adjustable Mobile Shower Commode is a versatile and convenient toileting and bathing solution designed for individuals with mobility limitations. This innovative commode features a height-adjustable design, making it suitable for users of different heights and needs.

The KCare Mobile Shower Commode is constructed from high-quality materials that ensure long-lasting use and easy maintenance. It has a sturdy frame that can support up to 200 kg and features adjustable legs that can be easily adjusted to accommodate different users' needs. The seat and backrest are padded and covered in a soft, comfortable material that provides additional comfort and support.

FEATURES

- 445mm seat width.
- Height adjustability of the frame.
- Flip up arms for unrestricted side transfer.
- Strap back rest for comfort and easy access.
- Durable and long lasting with a stainless steel frame.
- Swing away footrests.
- 4 x 125mm locking swivel castors.
- Conforms to AS/NZS 3973:2009 for safety and peace of mind.
- **Seat, pan, and pan carrier not included - sold separately**

SPECIFICATIONS

Max User Weight	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Frame Height Front Range	Floor to Seat Frame Height Rear Range
200kg	445mm	430mm	400mm	520mm	905mm	575 - 475mm	552 - 452mm