



KCARE TRANSIT MOBILE SHOWER COMMODE - HEIGHT ADJUSTABLE

\$1,598.00 Original price was: \$1,598.00.<u>\$1,438.00</u>Current price is: \$1,438.00.



SKU: KA114S-HA

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The KCare Mobile Height Adjustable Mobile Shower Commode is a versatile and convenient toileting and bathing solution designed for individuals with mobility limitations. This innovative commode features a height-adjustable design, making it suitable for users of different heights and needs.

The KCare Mobile Shower Commode is constructed from high-quality materials that ensure long-lasting use and easy maintenance. It has a sturdy frame that can support up to 200 kg and features adjustable legs that can be easily adjusted to accommodate different users' needs. The seat and backrest are padded and covered in a soft, comfortable material that provides additional comfort and support.

FEATURES

- 445mm seat width.
- Height adjustability of the frame.
- Flip up arms for unrestricted side transfer.
- Strap back rest for comfort and easy access.
- Durable and long lasting with a stainless steel frame.
- Swing away footrests.
- 4 x 125mm locking swivel castors.
- Conforms to AS/NZS 3973:2009 for safety and peace of mind.
- Seat, pan, and pan carrier not included sold separately



SPECIFICATIONS

Max User	Seat	Seat	Backrest	Overall	Overall	Frame Height	Frame Height
Weight	Width	Depth	Height	Width	Depth	Front Range	
200kg	445mm	430mm	400mm	520mm	905mm	575 - 475mm	552 - 452mm