




## KCARE MURRAY BRIDGE HIGH BACK CHAIR

**\$632.00 - \$775.00** Price range: \$632.00 through \$775.00



**SKU:** 18002HAB-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of KCare Murray Bridge High Back Chair

The KCare Murray Bridge High Back Chair is designed to provide superior comfort and support, especially for elderly users or those with mobility impairments. This chair is perfect for residential care facilities, hospitals, or home use, offering a high level of comfort and support with a stylish appearance.

### Features

- **Ultra Lightweight Design:** The chair features an ultra-lightweight aluminium frame with a timber look, combining strength with aesthetic appeal.
- **Sculptured Backrest:** Provides full lumbar support, enhancing comfort and aiding in proper posture alignment.
- **Removable Backrest:** For easy transport and adaptability in different settings.
- **Adjustable Height and Depth:** Offers customisation to fit various user needs, enhancing comfort and accessibility.
- **Variety of Seat Widths:** Available in 45 cm, 52 cm, and 60 cm seat widths, accommodating different body types.

### Specifications

- **Dimensions:** Comes in 45cm, 52cm and 60cm seat widths.
- **Colour Options:** Dot Coffee, Blue, Dot Forrest, Dot Ruby, Fawn Vinyl

Min Overall Depth

730 mm

Max Overall Depth	830 mm
Min Overall Height	1120 mm
Max Overall Height	1240 mm
Seat Width	600 mm
Min Seat Depth	430 mm
Max Seat Depth	530 mm
Min Seat Height	440 mm
Max Seat Height	560 mm
Backrest Height	720 mm
Product Weight	10.4 kg
Max User Weight	250 kg

The KCare Murray Bridge High Back Chair is an ideal solution for those who require extra support but also desire a piece of furniture that is both functional and stylish. Its adjustable features and comfortable design make it a valuable addition to any care setting or home environment.