


KCARE PADDED OPEN FRONT PU COMMODE SEAT - 445MM



SKU: KA11OPU

PRODUCT DESCRIPTION

FEATURES

- Polyurethane injection moulded seat.
- Includes pan carrier runners.
- Cost effective, comfortable economy seat.
- Clip-on design allows quick removal for cleaning and servicing.
- Suits 445mm wide mobile shower commodes only.
- Open and closed front seats available.

SPECIFICATIONS

Overall Width	440 mm
Overall Depth	430 mm
Overall Height	40 mm
Aperture Width	200 mm
Product Weight	3 kg
Brand	K Care

View the [KCare Tilt In Space Shower Commode](#) or our entire range of [shower commodes](#).

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au