




KCARE OVER TOILET FRAME 3 IN 1



SKU: KA450HDWA

PRODUCT DESCRIPTION

Purpose & Benefits of the Kcare Over Toilet Frame 3 in 1

The K Care 3-in-1 Width Adjustable Over Toilet Frame is a versatile aid designed to assist individuals who experience difficulty using standard toilets. Its multifunctional design allows it to serve as an over toilet aid, a freestanding commode, or a shower chair, providing flexibility to meet various user needs. The frame's width adjustability ensures a comfortable and secure fit for different body types and bathroom configurations.

Features

- **Width Adjustable Frame:** The frame can be adjusted to accommodate various user sizes and toilet dimensions, enhancing comfort and accessibility.
- **Height Adjustable Legs:** Equipped with height-adjustable legs to ensure the optimal seating height for each user, promoting ease of use and safety.
- **Durable Construction:** Constructed from white powder-coated 32 mm aluminium tubing, offering a lightweight yet robust structure that supports substantial weight while resisting corrosion.
- **Padded Seat with Open Front:** Features a removable padded seat with an open front design, providing comfort during use and facilitating easy cleaning and hygiene maintenance.
- **Removable Pan:** Includes a pan for use as a commode, which can be easily removed for cleaning or when the frame is positioned over a toilet.
- **Non-Slip Rubber Tips:** Fitted with non-slip rubber tips on each leg to prevent movement and enhance stability on various floor surfaces

Specifications

Min Overall Width	600 mm
Max Overall Width	750 mm
Min Overall Height	677 mm
Max Overall Height	777 mm
Seat Width	460 mm
Seat Depth	415 mm
Min Seat Height	500 mm
Max Seat Height	600 mm
Product Weight	4.6 kg
Max User Weight	300 kg

The Kcare Over Toilet Frame 3 in 1 is an essential aid for individuals requiring additional support, offering a safe and adaptable solution for toileting and hygiene needs.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au