



## KCARE OVER TOILET FRAME WITH SEAT FLAP - WHITE



**SKU:** KA410Z

## PRODUCT DESCRIPTION

### Purpose & Benefits of the KCare Over Toilet Frame with Seat Flap

The KCare Over Toilet Frame with Seat Flap is designed to provide additional support and safety for individuals who require assistance while toileting. Its adjustable features and sturdy construction make it suitable for both institutional and home use.

### Features

- **Adjustable Height and Angle** – Allows customisation to suit individual users and accommodate uneven floors, enhancing comfort and stability.
- **Clip-On Seat and Removable Splash Guard** – Facilitates easy assembly, cleaning, and maintenance, promoting hygiene.
- **Anti-Slip Rubber Tips** – Equipped on each leg to provide stability and prevent slipping during use.
- **Durable Construction** – Made from white powder-coated 22mm steel tubing, ensuring long-lasting durability.
- **Compliance with Standards** – Independently tested to AS/NZS 3973:2009, ensuring adherence to safety and quality standards.

### Specifications

- **Overall Width:** 510 mm
- **Overall Depth:** 600 mm
- **Minimum Overall Height:** 680 mm
- **Maximum Overall Height:** 805 mm
- **Seat Width:** 450 mm

- **Minimum Seat Height:** 480 mm
- **Maximum Seat Height:** 605 mm
- **Product Weight:** 5 kg
- **Maximum User Weight:** 110 kg

The KCare Over Toilet Frame with Seat Flap offers a practical and reliable solution for individuals seeking enhanced support and safety in toileting activities.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)