



KCARE OVER TOILET FRAME WITH SEAT FLAP

\$194.00



SKU: KA410AAHD

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The KCare Over Toilet Frame is designed for users requiring additional seating support while toileting. It has a white powder coat and 22mm steel tube legs. This toilet frame is height adjustable and can be angled to suit the individual users on uneven floors. The seat clips on and the splash guard removes easily for quick assembly, cleaning or replacing.

FEATURES

- Designed for users preferring additional seating support while toileting.
- Larger 28mm tubing makes it easier to grip handles when transferring.
- Adjustable height and angle to suit individual users and uneven floors.
- Seat clips on for easy assembly, change and cleaning.
- Supplied with toilet flap with removable splash guard.
- Anti-slip rubber tips on each leg for stability.
- Durable and long lasting aluminium frame.
- Independently tested to AS/NZS 3973:2009.
- Suitable for institutional and domestic use.

SPECIFICATIONS

Max User Weight	200 kg
Overall Width	530 mm
Min Overall Depth	570 mm
Max Overall Depth	605 mm

Seat Width	380 mm
Seat Depth	440 mm
Min Floor To Seat Height	470 mm
Max Floor To Seat Height	595 mm
Brand	K Care
Product Weight	2.8 kg