



KCARE OVER TOILET FRAME WITH SEAT FLAP



SKU: KA410AAHD

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The KCare Over Toilet Frame is designed for users requiring additional seating support while toileting. It has a white powder coat and 22mm steel tube legs. This toilet frame is height adjustable and can be angled to suit the individual users on uneven floors. The seat clips on and the splash guard removes easily for quick assembly, cleaning or replacing.

FEATURES

- Designed for users preferring additional seating support while toileting.
- Larger 28mm tubing makes it easier to grip handles when transferring.
- Adjustable height and angle to suit individual users and uneven floors.
- Seat clips on for easy assembly, change and cleaning.
- Supplied with toilet flap with removable splash guard.
- Anti-slip rubber tips on each leg for stability.
- Durable and long lasting aluminium frame.
- Independently tested to AS/NZS 3973:2009.
- Suitable for institutional and domestic use.

SPECIFICATIONS

Max User Weight	200 kg
Overall Width	530 mm
Min Overall Depth	570 mm
Max Overall Depth	605 mm

Seat Width	380 mm
Seat Depth	440 mm
Min Floor To Seat Height	470 mm
Max Floor To Seat Height	595 mm
Brand	K Care
Product Weight	2.8 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au