



## KCARE PADDED SHOWER STOOL WITH ARMS AND PADDED SEAT

**\$261.00**



**SKU: KA22210**

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

Designed for use in the shower, this KCare Padded Shower Stool has a padded backrest and seat. The padded seat is removable for cleaning. The seat height and angle are adjustable. Arms are provided for support while sitting or transferring. Rubber tips on each leg provide stability. The shower stool has an epoxy powder coated zinc treated mild steel, white frame and water proof polyurethane seat and back rest.

## FEATURES

- Provides a rest for the user while showering and is ideal in a smaller shower area
- Adjustable seat height and angle facilitates a comfortable seating position
- The clip on seat allows for easy assembly, cleaning and replacement
- Padded Seat option also available
- This shower stool also has the addition of a fixed padded backrest
- The backrest is permanently fixed to the stool frame
- Zinc treated mild steel frame with a powder coated finish

## SPECIFICATIONS

<b>Max Seat Height</b>	640mm
<b>Max User Weight</b>	110kg
<b>Min Seat Height</b>	520mm
<b>Overall Depth</b>	455mm
<b>Overall Width</b>	510mm

<b>Product Weight</b>	4kg
<b>Seat Depth</b>	310mm
<b>Seat Width</b>	360mm