



## KCARE PADDED SHOWER STOOL WITH ARMS AND PADDED SEAT



**SKU:** KA22210

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Designed for use in the shower, this KCare Padded Shower Stool has a padded backrest and seat. The padded seat is removable for cleaning. The seat height and angle are adjustable. Arms are provided for support while sitting or transferring. Rubber tips on each leg provide stability. The shower stool has an epoxy powder coated zinc treated mild steel, white frame and water proof polyurethane seat and back rest.

### FEATURES

- Provides a rest for the user while showering and is ideal in a smaller shower area
- Adjustable seat height and angle facilitates a comfortable seating position
- The clip on seat allows for easy assembly, cleaning and replacement
- Padded Seat option also available
- This shower stool also has the addition of a fixed padded backrest
- The backrest is permanently fixed to the stool frame
- Zinc treated mild steel frame with a powder coated finish

### SPECIFICATIONS

<b>Max Seat Height</b>	640mm
<b>Max User Weight</b>	110kg
<b>Min Seat Height</b>	520mm
<b>Overall Depth</b>	455mm
<b>Overall Width</b>	510mm

<b>Product Weight</b>	4kg
<b>Seat Depth</b>	310mm
<b>Seat Width</b>	360mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)