



KCARE PREMIUM KNEE WALKER



SKU: KA337

PRODUCT DESCRIPTION

Purpose & Benefits of the Kcare Premium Knee Walker

The Kcare Premium Knee Walker is a mobility solution designed to assist individuals recovering from foot, ankle, or lower leg injuries. It provides a comfortable, stable, and efficient way to maintain mobility without the need for crutches. With its adjustable design and ergonomic knee pad, the knee walker ensures a smooth and comfortable ride, allowing users to navigate both indoor and outdoor spaces with ease. Ideal for short-term rehabilitation, the Kcare Premium Knee Walker helps users maintain their independence while recovering.

Key Features

- **Ergonomic Knee Pad**: The contoured knee pad offers comfort and support, reducing pressure on the injured leg during movement, making it suitable for extended use.
- Adjustable Handlebar and Knee Pad Height: Easily adjust the handlebar and knee pad height to ensure a customised fit for users of various heights, promoting proper posture and comfort.
- **Durable, Lightweight Frame**: Constructed from high-quality materials, the lightweight frame provides robust support while remaining easy to manoeuvre.
- Large, All-Terrain Wheels: Equipped with large wheels that offer stability and smooth navigation over various surfaces, making it ideal for indoor and outdoor use.
- **Hand-Operated Brake System**: The built-in hand brakes ensure easy control and safe stopping, providing added security and peace of mind during movement.

Specifications

• Knee Pad Dimensions: 35 cm (L) x 18 cm (W)



Handlebar Height Range: Adjustable from 84 cm to 100 cm
Knee Pad Height Range: Adjustable from 48 cm to 58 cm

• Weight Capacity: 136 kg

• Frame Material: Lightweight aluminium frame

• Wheel Size: 20 cm diameter

• Cleaning: Wipe clean with a damp cloth

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 319903 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468

support@medimarnt.com.au