



## KCARE PREMIUM KNEE WALKER

**\$442.00**



**SKU: KA337**

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Kcare Premium Knee Walker

The Kcare Premium Knee Walker is a mobility solution designed to assist individuals recovering from foot, ankle, or lower leg injuries. It provides a comfortable, stable, and efficient way to maintain mobility without the need for crutches. With its adjustable design and ergonomic knee pad, the knee walker ensures a smooth and comfortable ride, allowing users to navigate both indoor and outdoor spaces with ease. Ideal for short-term rehabilitation, the Kcare Premium Knee Walker helps users maintain their independence while recovering.

### Key Features

- **Ergonomic Knee Pad:** The contoured knee pad offers comfort and support, reducing pressure on the injured leg during movement, making it suitable for extended use.
- **Adjustable Handlebar and Knee Pad Height:** Easily adjust the handlebar and knee pad height to ensure a customised fit for users of various heights, promoting proper posture and comfort.
- **Durable, Lightweight Frame:** Constructed from high-quality materials, the lightweight frame provides robust support while remaining easy to manoeuvre.
- **Large, All-Terrain Wheels:** Equipped with large wheels that offer stability and smooth navigation over various surfaces, making it ideal for indoor and outdoor use.
- **Hand-Operated Brake System:** The built-in hand brakes ensure easy control and safe stopping, providing added security and peace of mind during movement.

### Specifications

- **Knee Pad Dimensions:** 35 cm (L) x 18 cm (W)

- **Handlebar Height Range:** Adjustable from 84 cm to 100 cm
- **Knee Pad Height Range:** Adjustable from 48 cm to 58 cm
- **Weight Capacity:** 136 kg
- **Frame Material:** Lightweight aluminium frame
- **Wheel Size:** 20 cm diameter
- **Cleaning:** Wipe clean with a damp cloth