



KCARE PREMIUM PADDED SHOWER STOOL 520MM

~~\$304.00~~ Original price was: \$304.00. \$273.60 Current price is: \$273.60.



SKU: 12325A52

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **KCare Premium Padded Shower Stool (520mm)** is designed to provide comfort, safety, and stability for individuals who experience difficulty standing while showering. Its soft PU padded seat offers superior comfort, while the durable 32mm aluminium frame ensures long-lasting strength and resistance to corrosion. Ideal for everyday use, this stool promotes independence and confidence in the bathroom.

FEATURES

- Assists users who have difficulty showering independently
- Soft PU padded seat provides enhanced comfort and support
- Height-adjustable legs for personalised fit and optimal stability
- Constructed from 32mm powder-coated aluminium tubing for strength and durability
- Lightweight, rust-resistant frame suitable for wet environments
- Available in 450mm, 520mm, and 600mm frame widths

SPECIFICATIONS

- **Min Overall Height:** 670 mm
- **Max Overall Height:** 810 mm
- **Seat Width:** 520 mm
- **Seat Depth:** 300 mm
- **Seat Height Range:** 480 – 640 mm
- **Max User Weight:** 250 kg

- **Frame Material:** Powder-coated 32mm aluminium tubing
- **Seat Material:** PU padded
- **Brand:** KCare