



## KCARE REHAB PAEDIATRIC SELF PROPELLED SHOWER COMMODE



**SKU:** KA120SRC-REHAB-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Kcare Rehab Shower Commode Self Propelled - Paediatric, a cutting-edge solution designed to aid the mobility and comfort of children who require support in daily showering and toileting routines. Engineered with versatility and longevity in mind, this product combines an array of features to adapt to the unique needs of each user.

### FEATURES

- **Seat Widths:** 360 mm and 320 mm options for a tailored fit.
- **Lateral Thoracic Support:** Adjustable height mounting plates with 240 mm range.
- **Height Adjustability:** Adapts to the growing child.
- **Knee and Hip Guides:** Standard mounting plates for alignment.
- **Unrestricted Side Transfer:** Flip-up arms and flush rear wheels.
- **Comfortable Backrest:** Adjustable tension for support.
- **Flexible Seat Frame:** Left and right side opening options (not standard).
- **Seat Depth Adjustment:** Customisable for comfort.
- **Durability:** Stainless steel frame for longevity.
- **Footrests:** Swing-away design for convenience.
- **Mobility Support:** 2 x Front 125 mm locking swivel castors.
- **[Seat](#), [pan](#), and [pan carrier](#) not included - sold separately**

## SPECIFICATIONS

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Floor to		Max User Weight
						Seat Height Front Range	Seat Height Rear Range	
Child Self-Propelled Rehab Commode KA120SRK-REHAB	360mm	320mm	360mm	630mm	900mm	550mm 475mm	527mm 452mm	120kg
Paediatric Self-Propelled Rehab Commode KA120SRC-REHAB	330mm	320mm	360mm	630mm	900mm	550mm 475mm	527mm 452mm	120kg

## BROCHURE

[KCare\\_Rehab\\_Shower\\_Commode\\_Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)