




KCARE REHAB SELF PROPELLED HEIGHT ADJUSTABLE SHOWER COMMODE



SKU: KA120SR-REHAB-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Kcare Rehab Shower Commode Self-Propelled - a leading-edge solution designed for ultimate comfort and accessibility. This robust shower commode combines functionality, adjustability, and durability, creating a user-centred experience for individuals with special needs.

FEATURES

- **Seat Widths:** Available in 400mm, 445mm, and 500mm widths for optimal fit.
- **Lateral Support:** Adjustable thoracic support with a 240mm range.
- **Height Adjustability:** Both frame and flip-up arms can be adjusted for personalised comfort.
- **Knee and Hip Guides:** Standard on all models with a 40mm horizontal range.
- **Adjustable Backrest:** Tailor the tension for your comfort needs.
- **Versatile Seat Frame:** Standard openings on both sides for added flexibility.
- **Seat Depth Adjustment:** Easy customisation for enhanced sitting experience.
- **Sturdy Construction:** Durable stainless steel frame with long-lasting build.
- **Swing Away Footrests:** For added comfort in positioning.
- **Self-Propelling Design:** Large rear wheels facilitate unrestricted side transfer.
- **Seat, pan, and pan carrier not included - sold separately**

SPECIFICATIONS

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Height Front Range	Floor to Seat Height Rear Range	Max User Weight
Junior Self Propelled Rehab Commode KA120SRJ-REHAB	400mm	435mm	360mm	585mm	1050mm	-	-	200kg
Adult Self Propelled Rehab Commode KA120SR-REHAB	445mm	435mm	360mm	630mm	1050mm	-	-	200kg
Adult Self Propelled Rehab Commode KA120SR50-REHAB	500mm	435mm	360mm	685mm	1050mm	-	-	200kg

BROCHURE

[KCare_Rehab_Shower_Commode_Brochure](#)