



KCARE REHAB SELF PROPELLED HEIGHT ADJUSTABLE SHOWER COMMODE



SKU: KA120SR-REHAB-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Kcare Rehab Shower Commode Self-Propelled - a leading-edge solution designed for ultimate comfort and accessibility. This robust shower commode combines functionality, adjustability, and durability, creating a user-centred experience for individuals with special needs.

FEATURES

- Seat Widths: Available in 400mm, 445mm, and 500mm widths for optimal fit.
- Lateral Support: Adjustable thoracic support with a 240mm range.
- Height Adjustability: Both frame and flip-up arms can be adjusted for personalised comfort.
- Knee and Hip Guides: Standard on all models with a 40mm horizontal range.
- Adjustable Backrest: Tailor the tension for your comfort needs.
- Versatile Seat Frame: Standard openings on both sides for added flexibility.
- Seat Depth Adjustment: Easy customisation for enhanced sitting experience.
- Sturdy Construction: Durable stainless steel frame with long-lasting build.
- **Swing Away Footrests**: For added comfort in positioning.
- Self-Propelling Design: Large rear wheels facilitate unrestricted side transfer.
- Seat, pan, and pan carrier not included sold separately

SPECIFICATIONS



Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth		Floor to Seat Height Rear Range	Max User Weight
Junior Self Propelled		435mm				550mm	527mm	
Rehab Commode	400mm	-	360mm	585mm	1050mm	-	-	200kg
KA120SRJ-REHAB		465mm				475mm	452mm	
Adult Self Propelled		435mm				550mm	527mm	
Rehab Commode	445mm	-	360mm	630mm	1050mm	-	_	200kg
KA120SR-REHAB		465mm				475mm	452mm	
Adult Self Propelled		435mm				550mm	527mm	
Rehab Commode	500mm	-	360mm	685mm	1050mm	-	_	200kg
KA120SR50-REHAB		465mm				475mm	452mm	

BROCHURE

KCare_Rehab_Shower_Commode_Brochure

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom	Elizabeth Vale Showroom	Reynella Showroom
343 Brighton Rd, North	Playford Health Hub, 44 John Rice	Shop 1/141 Sherriffs Rd,
Brighton SA 5048	Ave, Elizabeth Vale SA 5112	Reynella SA 5161
08 7079 8773	08 7079 8777	08 7079 8771
brighton@medimart.com.au	elizabethvale@medimart.com.au	reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au



NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au