



KCARE REHAB TILT IN SPACE SHOWER COMMODE



SKU: KA11507-REHAB-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Kcare Rehab Tilt In Space Shower Commode with Swingaway Footrest offers superior functionality and adjustability to meet the needs of patients with varying sizes and conditions. Available in two seat widths - 445mm and 500mm - this commode is designed to provide an optimal fit for individual comfort.

Featuring lateral thoracic support adjustable height mounting plates standard on all models, it allows a 240mm range (20mm increments) for precise customisation. Knee guide mounting plates, also standard on all models, provide a 40mm horizontal mounting range (20mm increments), adding to its adaptability.

The height adjustability of the frame allows further personalisation, while height-adjustable flip-up arms enable unrestricted side transfers. Comfort is enhanced by adjustable tension backrest upholstery and the inclusion of standard left and right side openings on the seat frame.

Built to last with a durable stainless steel frame, the K Care Tilt in Space Rehab Shower Commode is both robust and long-lasting. Swing-away footrests with padded calf supports provide added comfort and ease of use, and 4 x 125mm locking swivel castors ensure mobility and stability.

FEATURES

- 445mm and 500mm seat widths
- Lateral thoracic support adjustable height mounting plates are standard on all models. 240mm range (20mm increments)
- Height adjustability of the frame

- Knee guide mounting plates are standard on all models. 40mm horizontal mounting range (20mm increments)
- Height adjustable flip-up arms for unrestricted side transfer
- Adjustable tension backrest upholstery
- Standard left and right side openings on the seat frame
- Durable and long-lasting with a stainless steel frame
- Swing-away footrests with padded calf support
- 4 x 125mm locking swivel castors
- **Seat, pan, and pan carrier not included - sold separately**

SPECIFICATIONS

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Max User Weight
KA11507-REHAB 445mm		435mm - 465mm	855mm	580mm	580mm	200kg
KA11512-REHAB 500mm		435mm - 465mm	855mm	595mm	580mm	200kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au