




KCARE REHAB PAEDIATRIC TRANSIT SHOWER COMMODE



SKU: KA114SC-REHAB-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Kcare Rehab Shower Commode Transit - Paediatric, a specialised and innovative solution designed to enhance comfort and support for young users. This state-of-the-art shower commode has been meticulously crafted to meet various needs, providing a personalised experience that combines functionality with durability.

Available in both 320mm and 360mm seat widths, it ensures a snug and comfortable fit. The lateral thoracic support features adjustable height mounting plates standard on all models, with a 240mm range in 20mm increments, providing tailored support for the upper body.

With the height adjustability of the frame, it's adaptable to various body sizes and shapes. Knee and hip guide mounting plates offer a 40mm horizontal mounting range in 20mm increments, allowing precise alignment and support for the legs. The height-adjustable flip-up arms allow easy and unrestricted side transfer, adding convenience for both users and caregivers.

FEATURES

- **360mm and 320mm Seat Widths:** Available for a comfortable fit.
- **Lateral Thoracic Support:** Adjustable height mounting plates are standard on all models, with a 240mm range (20mm increments) for tailored upper body support.
- **Height Adjustability of the Frame:** Adaptable to various body sizes and shapes.
- **Knee and Hip Guide Mounting Plates:** Standard on all models, offering a 40mm horizontal mounting range (20mm increments) for precise leg alignment.
- **Height Adjustable Flip-Up Arms:** Allows for unrestricted side transfer, enhancing user

convenience.

- **Adjustable Tension Backrest Upholstery:** Ensures comfort by adapting to individual body contours.
- **Option of Left and Right Side Opening on Seat Frame:** Provides flexibility in positioning (not standard on this model).
- **Seat Depth Adjustment:** Offers additional customisation for comfort.
- **Durable and Long-Lasting Construction:** Built with a stainless steel frame for long-lasting use.
- **Swing Away Footrests:** Enhances ease of access and use.
- **4 x 125mm Locking Swivel Castors:** Ensures easy mobility and stable locking as required.
- **[Seat](#), [pan](#), and [pan carrier](#) not included - sold separately**

These features combine to make the Kcare Rehab Shower Commode Transit - Paediatric a versatile and user-friendly solution for those in need of specialised care.

SPECIFICATIONS

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Height Front Range	Floor to Seat Height Rear Range	Max User Weight
Child Transit Rehab Commode KA114SK-REHAB	360mm	320mm	360mm	485mm	750mm	575mm - 475mm	552mm - 452mm	120kg
Paediatric Transit Rehab Commode KA114SC-REHAB	320mm	320mm	360mm	485mm	750mm	575mm - 475mm	552mm - 452mm	120kg

BROCHURE

[KCare_Rehab_Shower_Commode_Brochure](#)