




## KCARE REHAB PAEDIATRIC TRANSIT SHOWER COMMODE



**SKU:** KA114SC-REHAB-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Kcare Rehab Shower Commode Transit - Paediatric, a specialised and innovative solution designed to enhance comfort and support for young users. This state-of-the-art shower commode has been meticulously crafted to meet various needs, providing a personalised experience that combines functionality with durability.

Available in both 320mm and 360mm seat widths, it ensures a snug and comfortable fit. The lateral thoracic support features adjustable height mounting plates standard on all models, with a 240mm range in 20mm increments, providing tailored support for the upper body.

With the height adjustability of the frame, it's adaptable to various body sizes and shapes. Knee and hip guide mounting plates offer a 40mm horizontal mounting range in 20mm increments, allowing precise alignment and support for the legs. The height-adjustable flip-up arms allow easy and unrestricted side transfer, adding convenience for both users and caregivers.

### FEATURES

- **360mm and 320mm Seat Widths:** Available for a comfortable fit.
- **Lateral Thoracic Support:** Adjustable height mounting plates are standard on all models, with a 240mm range (20mm increments) for tailored upper body support.
- **Height Adjustability of the Frame:** Adaptable to various body sizes and shapes.
- **Knee and Hip Guide Mounting Plates:** Standard on all models, offering a 40mm horizontal mounting range (20mm increments) for precise leg alignment.
- **Height Adjustable Flip-Up Arms:** Allows for unrestricted side transfer, enhancing user

convenience.

- **Adjustable Tension Backrest Upholstery:** Ensures comfort by adapting to individual body contours.
- **Option of Left and Right Side Opening on Seat Frame:** Provides flexibility in positioning (not standard on this model).
- **Seat Depth Adjustment:** Offers additional customisation for comfort.
- **Durable and Long-Lasting Construction:** Built with a stainless steel frame for long-lasting use.
- **Swing Away Footrests:** Enhances ease of access and use.
- **4 x 125mm Locking Swivel Castors:** Ensures easy mobility and stable locking as required.
- **Seat, pan, and pan carrier not included - sold separately**

These features combine to make the Kcare Rehab Shower Commode Transit - Paediatric a versatile and user-friendly solution for those in need of specialised care.

## SPECIFICATIONS

| Description                                    | Seat Width | Seat Depth | Backrest Height | Overall Width | Overall Depth | Floor to Seat Height Front Range | Floor to Seat Height Rear Range | Max User Weight |
|--|------------|------------|-----------------|---------------|---------------|----------------------------------|---------------------------------|-----------------|
| Child Transit Rehab Commode KA114SK-REHAB      | 360mm      | 320mm      | 360mm           | 485mm         | 750mm         | 575mm - 475mm                    | 552mm - 452mm                   | 120kg           |
| Paediatric Transit Rehab Commode KA114SC-REHAB | 320mm      | 320mm      | 360mm           | 485mm         | 750mm         | 575mm - 475mm                    | 552mm - 452mm                   | 120kg           |

## BROCHURE

[KCare\\_Rehab\\_Shower\\_Commode\\_Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

**Brighton Showroom**

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

**Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

**Reynella Showroom**

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

**VIC****Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

**Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

**NT****Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)