



KCARE REHAB TRANSIT SHOWER COMMODE



SKU: KA114S-REHAB-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Kcare Rehab Transit Shower Commode, an innovative commode solution tailored for comfort and support. Available in 400mm, 445mm, and 500mm seat widths, it ensures a perfect fit. Lateral thoracic support, adjustable height mounting plates, and knee and hip guide mounting plates provide personalised alignment and comfort.

The frame's height adjustability, flip-up arms, adjustable tension backrest upholstery, and seat depth adjustment offer maximum flexibility. Standard left and right side openings on the seat frame enhance practicality, while the durable stainless steel frame ensures long-lasting use. Swing-away footrests and 4 x 125mm locking swivel castors add to the convenience and mobility.

FEATURES

- **400mm, 445mm, and 500mm seat widths:** Choose the right size for comfort.
- **Lateral thoracic support adjustable height mounting plates:** Standard on all models, with a 240mm range (20mm increments).
- **Height adjustability of the frame:** Customise to the user's height.
- **Knee and Hip guide mounting plates:** Standard on all models, with a 40mm horizontal mounting range (20mm increments).
- **Height adjustable flip-up arms:** Allow unrestricted side transfer.
- **Adjustable tension backrest upholstery:** Ensures comfort in every position.
- **Standard left and right side opening on seat frame:** Enhances practicality.
- **Seat depth adjustment:** Accommodates various user requirements.
- **Durable and long-lasting with a stainless steel frame:** For extended use.
- **Swing-away footrests:** Provide added comfort and accessibility.

- **4 x 125mm locking swivel castors:** Ensure stability and easy movement.
- **[Seat](#), [pan](#), and [pan carrier](#) not included - sold separately**

SPECIFICATIONS

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Height (Front Range)	Floor to Seat Height (Rear Range)	Max User Weight
Junior Transit Rehab Commode KA114SJ-REHAB	400mm	435mm	360mm	500mm	905mm	575mm - 475mm	552mm - 452mm	200kg
Adult Transit Rehab Commode KA114S-REHAB	445mm	435mm	360mm	540mm	905mm	575mm - 475mm	552mm - 452mm	200kg
Adult Transit Rehab Commode KA114S50-REHAB	500mm	435mm	360mm	595mm	905mm	575mm - 475mm	552mm - 452mm	200kg

BROCHURE

[KCare_Rehab_Shower_Commode_Brochure](#)