



KCARE REHAB TRANSIT SHOWER COMMODE



SKU: KA114S-REHAB-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Kcare Rehab Transit Shower Commode, an innovative commode solution tailored for comfort and support. Available in 400mm, 445mm, and 500mm seat widths, it ensures a perfect fit. Lateral thoracic support, adjustable height mounting plates, and knee and hip guide mounting plates provide personalised alignment and comfort.

The frame's height adjustability, flip-up arms, adjustable tension backrest upholstery, and seat depth adjustment offer maximum flexibility. Standard left and right side openings on the seat frame enhance practicality, while the durable stainless steel frame ensures long-lasting use. Swing-away footrests and 4 x 125mm locking swivel castors add to the convenience and mobility.

FEATURES

- **400mm, 445mm, and 500mm seat widths:** Choose the right size for comfort.
- **Lateral thoracic support adjustable height mounting plates:** Standard on all models, with a 240mm range (20mm increments).
- **Height adjustability of the frame:** Customise to the user's height.
- **Knee and Hip guide mounting plates:** Standard on all models, with a 40mm horizontal mounting range (20mm increments).
- **Height adjustable flip-up arms:** Allow unrestricted side transfer.
- **Adjustable tension backrest upholstery:** Ensures comfort in every position.
- **Standard left and right side opening on seat frame:** Enhances practicality.
- **Seat depth adjustment:** Accommodates various user requirements.
- **Durable and long-lasting with a stainless steel frame:** For extended use.
- **Swing-away footrests:** Provide added comfort and accessibility.

- **4 x 125mm locking swivel castors:** Ensure stability and easy movement.
- **Seat, pan, and pan carrier not included - sold separately**

SPECIFICATIONS

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Height (Front Range)	Floor to Seat Height (Rear Range)	Max User Weight
Junior Transit Rehab Commode KA114SJ-REHAB	400mm	435mm - 465mm	360mm	500mm	905mm	575mm - 475mm	552mm - 452mm	200kg
Adult Transit Rehab Commode KA114S-REHAB	445mm	435mm - 465mm	360mm	540mm	905mm	575mm - 475mm	552mm - 452mm	200kg
Adult Transit Rehab Commode KA114S50-REHAB	500mm	435mm - 465mm	360mm	595mm	905mm	575mm - 475mm	552mm - 452mm	200kg

BROCHURE

[KCare_Rehab_Shower_Commode_Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au