



## KCARE SELF PROPELLED HEIGHT ADJUSTABLE MOBILE SHOWER COMMODE



**SKU:** KA120SR-HA

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The KCare Self Propelled Height Adjustable Mobile Shower Commode is a versatile and practical solution for individuals with mobility challenges who require assistance with personal hygiene and toileting. This commode is designed to offer the user maximum comfort and convenience, while maintaining their independence.

This mobile shower commode features height-adjustable legs that can be easily adjusted to suit the user's height, ensuring a comfortable and safe transfer to and from the commode. The self-propelled function allows the user to move around freely and safely, without requiring assistance from a caregiver. The padded seat and backrest provide additional comfort during extended periods of use.

The KCare Self Propelled Height Adjustable Mobile Shower Commode is also equipped with a removable bucket, making it easy to clean and maintain. The chair is constructed from durable materials and is built to withstand frequent use and movement.

**Please note: Padded seat is not included**

### FEATURES

- Height adjustable for user comfort.
- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- 3 strap backrest for comfort and easy access.
- Integrated push-handle for manoeuvrability.
- Swing-up arms for unrestricted side transfer.
- Swingaway and removable height adjustable footrest.

- 2x 125mm total locking castors for stability.
- 2x 22" self-propel wheels inc. push to lock brakes for ease of movement.
- Designed to be used with a K Care patented removable seat.
- Push rims available as optional extra.
- **Seat, pan, and pan carrier not included - sold separately**

## SPECIFICATIONS

Overall Width	730 mm
Overall Depth	905 mm
Seat Width	445 mm
Seat Depth	430 mm
Min Floor To Seat Height Front	475 mm
Max Floor To Seat Height Front	550 mm
Min Floor To Seat Height Rear	452 mm
Max Floor To Seat Height Rear	527 mm
Armrest Type	Loop
Armrest Height	285 mm
Backrest Height	400 mm
Footrest Type	Swing Away
Product Weight	12 kg
Max User Weight	200 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

**Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

**NT**

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)