



# KCARE SELF PROPELLED MOBILE SHOWER COMMODE WITH SWINGAWAY FOOTREST

\$2,176.00



SKU: KA120SR-H

#### PRODUCT DESCRIPTION

### **PURPOSE & BENEFITS**

The Kcare Self Propel Mobile Shower Commode with Swingaway Footrest is a high-quality solution for individuals with mobility limitations who require assistance with toileting and personal hygiene. This commode is designed to offer maximum comfort and convenience while maintaining user independence.

The swingaway footrests provide additional support and stability during transfers and can be easily removed for cleaning or when not required. The chair also features a self-propelled function that allows the user to move around freely and safely without assistance from a caregiver.

#### Please note: Padded seat is not included

#### **FEATURES**

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- 3 strap backrest for comfort and easy access.
- Integrated push-handle for manoeuvrability.
- Swing-up arms for unrestricted side transfer.
- Swingaway and removable height adjustable footrest.
- 2x 125mm total locking castors for stability.
- 2x 22" self-propel wheels inc. push to lock brakes for ease of movement.
- Designed to be used with a Kcare patented removable seat.
- Available in 320, 400, 445, 500, 550, 600, and 650mm seat widths to accommodate various body types.
- Push rims available as optional extra.



• Seat, pan, and pan carrier not included - sold separately

## **SPECIFICATIONS**

					Floor					Max
Seat	Overall	Overall	Overall	Seat	То	Armrest	Backrest	Footrest	<b>Product</b>	User
Width	Width	Depth	Height	Depth		Type	Height	Type	Weight	Weight
					Height					
320	630	1050	1005	320	550	Loop	400 mm	Swing	16 kg	120 kg
mm	mm	mm	mm	mm	mm	СООР	400 111111	Away	10 kg	120 Kg
400	630	1050	1005	430	550	_oop 400 mm	Swing	16 kg	200 kg	
mm	mm	mm	mm	mm	mm	СООР	400 11111	Away	10 kg	200 Kg
445	630	1050	1005	430	550	Swingup	400 mm	Swing	15 kg	200 kg
mm	mm	mm	mm	mm	mm	Loop	400 111111	Away	13 Kg	200 kg
500	685	1045	1005	430	550	Swingup	400 mm	Swing	17 kg	200 kg
mm	mm	mm	mm	mm	mm	Loop	100 111111	Away	17 kg	200 kg
550	735	1045	1005	460	550	Padded	400 mm	Swing	19 kg	300 kg
mm	mm	mm	mm	mm	mm	Removable	100 111111	Away	13 Kg	Soo kg
600	785	1045	1005	470	550	Padded	400 mm	Swing	20 kg	300 kg
mm	mm	mm	mm	mm	mm	Removable	400 111111	Away	20 Kg	Joo kg
650	835	1045	1005	480	550	Padded	400 mm	Swing	20 kg	300 kg
mm	mm	mm	mm	mm	mm	Removable		Away	20 Ng	Joo Ng