



KCARE SELF PROPELLED MOBILE SHOWER COMMODE WITH SWINGAWAY FOOTREST

\$1,657.00 - \$2,681.25 Price range: \$1,657.00 through \$2,681.25



SKU: KA120SR-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Kcare Self Propel Mobile Shower Commode with Swingaway Footrest is a high-quality solution for individuals with mobility limitations who require assistance with toileting and personal hygiene. This commode is designed to offer maximum comfort and convenience while maintaining user independence.

The swingaway footrests provide additional support and stability during transfers and can be easily removed for cleaning or when not required. The chair also features a self-propelled function that allows the user to move around freely and safely without assistance from a caregiver.

Please note: Padded seat is not included

FEATURES

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- 3 strap backrest for comfort and easy access.
- Integrated push-handle for manoeuvrability.
- Swing-up arms for unrestricted side transfer.
- Swingaway and removable height adjustable footrest.
- 2x 125mm total locking castors for stability.
- 2x 22" self-propel wheels inc. push to lock brakes for ease of movement.
- Designed to be used with a Kcare patented removable seat.
- Available in 320, 400, 445, 500, 550, 600, and 650mm seat widths to accommodate various body types.
- Push rims available as optional extra.

- [Seat](#), [pan](#), and [pan carrier](#) not included - sold separately

SPECIFICATIONS

Seat Width	Overall Width	Overall Depth	Overall Height	Seat Depth	Floor To Seat Height	Armrest Type	Backrest Height	Footrest Type	Product Weight	Max User Weight
320 mm	630 mm	1050 mm	1005 mm	320 mm	550 mm	Loop	400 mm	Swing Away	16 kg	120 kg
400 mm	630 mm	1050 mm	1005 mm	430 mm	550 mm	Loop	400 mm	Swing Away	16 kg	200 kg
445 mm	630 mm	1050 mm	1005 mm	430 mm	550 mm	Swingup Loop	400 mm	Swing Away	15 kg	200 kg
500 mm	685 mm	1045 mm	1005 mm	430 mm	550 mm	Swingup Loop	400 mm	Swing Away	17 kg	200 kg
550 mm	735 mm	1045 mm	1005 mm	460 mm	550 mm	Padded Removable	400 mm	Swing Away	19 kg	300 kg
600 mm	785 mm	1045 mm	1005 mm	470 mm	550 mm	Padded Removable	400 mm	Swing Away	20 kg	300 kg
650 mm	835 mm	1045 mm	1005 mm	480 mm	550 mm	Padded Removable	400 mm	Swing Away	20 kg	300 kg