



KCARE SHOWER CHAIR WITH ARMS AND PLASTIC SEAT



SKU: KA220ZA

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Designed for users preferring additional seating support while showering. The KCare Shower Chair has a white powder coated 22mm steel tube, adjustable height and angle to suit individual users on uneven floors, and a contoured seat with inbuilt handle and water drainage slots. Clip on seat for easy assembly.

FEATURES

- Designed for users preferring additional seating support while showering.
- White powder coated 22mm steel tube.
- Adjustable height and angle to suit individual users and uneven floors.
- Contoured seat with inbuilt handle and water drainage slots.
- Seat clips on for easy assembly, change and cleaning.
- Anti-slip rubber tips on each leg for stability.
- Durable and long lasting aluminium frame.
- Independently tested to AS/NZS 3973:2009.
- Suitable for institutional and domestic use.

SPECIFICATIONS

Overall Width

505 mm

Overall Depth	450 mm
Seat Width	450 mm
Seat Depth	430 mm
Min Seat Height	420 mm
Max Seat Height	550 mm
Product Weight	4 kg
Max User Weight	110 kg
Brand	K Care

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au