




## KCARE SHOWER CHAIR WITH ARMS AND PLASTIC SEAT

**\$107.00**



**SKU:** KA220ZA

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Designed for users preferring additional seating support while showering. The KCare Shower Chair has a white powder coated 22mm steel tube, adjustable height and angle to suit individual users on uneven floors, and a contoured seat with inbuilt handle and water drainage slots. Clip on seat for easy assembly.

### FEATURES

- Designed for users preferring additional seating support while showering.
- White powder coated 22mm steel tube.
- Adjustable height and angle to suit individual users and uneven floors.
- Contoured seat with inbuilt handle and water drainage slots.
- Seat clips on for easy assembly, change and cleaning.
- Anti-slip rubber tips on each leg for stability.
- Durable and long lasting aluminium frame.
- Independently tested to AS/NZS 3973:2009.
- Suitable for institutional and domestic use.

### SPECIFICATIONS

Overall Width

505 mm

Overall Depth	450 mm
Seat Width	450 mm
Seat Depth	430 mm
Min Seat Height	420 mm
Max Seat Height	550 mm
Product Weight	4 kg
Max User Weight	110 kg
Brand	K Care