



## KCARE SHOWER STOOL WITH ARMS - HEIGHT ADJUSTABLE, CLIP-ON SEAT

**\$87.00**



**SKU: KA222ZA**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The KCare Shower Stool with Arms has a plastic seat that can easily be removed for cleaning/replacement and the height and angle can be adjusted. There are arms provide support to users while sitting or transferring. Rubber tips on each leg provide anti-slip stability. This shower stool is durable and is available with an epoxy powder coated zinc treated mild steel, white or aluminium frame and moulded PVC seat.

### FEATURES

- Designed to assist those who have difficulty showering.
- Height adjustable legs .
- Plastic seat for ease of cleaning .
- Powder coated white 22mm steel tubing .

### SPECIFICATIONS

Overall Width	510 mm
Overall Depth	380 mm
Product Weight	3 kg
Max User Weight	110 kg
Seat Width	340 mm
Seat Depth	290 mm
Min Seat Height	510 mm

Max Seat Height

630 mm